



Validated Windows 10

Windows 10 Pro	1607	1703	1803	1809	1903	1909	2004	20H2	21H1	21H2
LabSolutions LC/GC	5.111	5.92, 5.111	5.96 ~ 5.99, 5.111	5.96SP3 ~ 5.101, 5.111	5.98 ~ 5.106, 5.111	5.99 ~ 5.111	5.101 ~ 5.111	5.103 ~ 5.111	5.110 ~ 5.111	
LabSolutions LCMS	5.89, 5.114	5.91, 5.114	5.95 ~ 5.97, 5.114	5.114	5.98SP1 ~ 5.99SP2, 5.114	5.114	5.109 ~ 5.114	5.109 ~ 5.114	5.113 ~ 5.114	5.114
LabSolutions DB	6.111	6.82, 6.111	6.86 ~ 6.92, 6.111	6.86SP3 ~ 6.101, 6.111	6.89 ~ 6.106, 6.111	6.92 ~ 6.111	6.101 ~ 6.111	6.103 ~ 6.111	6.110 ~ 6.111	
LabSolutions DB LCMS	6.72, 6.114	6.81, 6.114	6.86 ~ 6.87SP1, 6.114	6.114	6.91, 6.114	6.93, 6.114	6.109 ~ 6.114	6.109 ~ 6.114	6.114	6.114
LabSolutions DB GCMS								6.107 ~ 6.116	6.116	6.116
LabSolutions CS	6.82SP1	6.82SP1	6.86 ~ 6.86SP1	6.86SP2 ~ 6.108	6.88 ~ 6.108	6.90SP2, 6.108	6.102 ~ 6.115	6.108 ~ 6.115	6.115	6.115
LabSolutions Manager DB *1			1.86 ~ 1.88		1.88		1.102 ~ 1.108	1.108		
LabSolutions Manager DB (MIDR) *1			1.88		1.88		1.108 ~ 1.115	1.108 ~ 1.115	1.115	1.115
LabSolutions Manager CS *1			1.86SP1	1.86SP2 ~ 1.90	1.88 ~ 1.90		1.102 ~ 1.115	1.108 ~ 1.115	1.115	1.115
LabSolutions LCMS CL	1.10	1.20								
LabSolutions CL										1.40

Windows 10 Enterprise					1903	1909	2004	20H2	21H1	21H2
LabSolutions LC/GC										
LabSolutions LCMS					5.114	5.114	5.114	5.114	5.114	5.114
LabSolutions DB					6.106, 6.111	6.106 ~ 6.111	6.106 ~ 6.111	6.106 ~ 6.111	6.110 ~ 6.111	
LabSolutions DB LCMS							6.109 ~ 6.114	6.109 ~ 6.114	6.114	6.114
LabSolutions DB GCMS								6.107 ~ 6.116	6.116	6.116
LabSolutions CS								6.108 ~ 6.115	6.115	6.115
LabSolutions Manager DB *1								1.108		
LabSolutions Manager DB (MIDR) *1								1.108 ~ 1.115	1.115	1.115
LabSolutions Manager CS *1								1.108 ~ 1.115	1.115	1.115

Windows 10 Enterprise LTSC		LTSC 2019 (1809)							LTSC 2021 (21H2)
LabSolutions LC/GC		5.106							
LabSolutions LCMS		5.113 ~ 5.114							5.114
LabSolutions DB		6.106, 6.111							
LabSolutions DB LCMS		6.109 ~ 6.114							6.114
LabSolutions DB GCMS		6.107 ~ 6.116							6.116
LabSolutions CS		6.86SP2 ~ 6.108							6.115
LabSolutions UV-Vis		1.12 ~ 1.13							
LabSolutions IR		2.30							
LabSolutions ICPMS		1.60 ~ 1.61							

--	--	--	--	--	--	--	--	--

Windows Edition	1803	1809	1903	1909	2004	20H2	21H1	21H2
UVProbe - Extension Program					2.71 - 1.04	2.71 - 1.04	2.71 - 1.04	2.71 - 1.04
LabSolutions UV-Vis - Analysis Sequence					1.11 - 1.12	1.12 - 1.12	1.12 ~ 1.13 - 1.12	1.12 ~ 1.13 - 1.12
UV Performance Validation					1.22	1.23	1.30	1.30
LabSolutions IR - Analysis Sequence - AIMsolution					2.27 - 1.12 - 1.25	2.27 - 1.12 - 1.25	2.30 - 1.12 - 1.30	2.30 - 1.12 - 1.30
LabSolutions RF - Extension Program - Analysis Sequence					1.17 - 1.04 - 1.12	1.17 - 1.04 - 1.12	1.17 - 1.04 - 1.12	1.17 - 1.04 - 1.12
LabSolutions ICPMS - Extension Program					1.50 - 1.50	1.50 ~ 1.61 - 1.50	1.60 ~ 1.61 - 1.60	1.60 ~ 1.61 - 1.60
WizAArd(AA-7800) - Agent Connection Kit					6.00 - 6.00	6.00 - 6.00	6.00 - 6.00	6.00 - 6.00
WizAArd(AA-7000) - Agent Connection Kit					5.24 / 5.25 - 5.10	5.24 / 5.25 - 5.10	5.24 / 5.25 - 5.10	5.24 / 5.25 - 5.10
ICPEsolution 9000 - MIDR Connection Kit					1.22 - 1.00	1.22 - 1.00	1.22 - 1.00	1.22 - 1.00
ICPEsolution 9800 - MIDR Connection Kit					2.22 - 1.00	2.22 - 1.00	2.22 - 1.00	2.22 - 1.00

TOC

Windows 10 Pro	1803	1809	1903	1909	2004	20H2	21H1	21H2
LabSolutions TOC		1.00	1.00	1.00	1.00	1.00	1.00	1.00
TOC-Control L	1.09	1.09	1.09	1.09	1.09	1.09		
TOC-Control V	2.60	2.60						

*1 UVやIRなどの各機種アプリと組み合わせて使用する場合、対応状況は各BUへご確認ください。
Please check with each BU for support of the software when you use each instrument specific software (e.g. UV-Vis, LabSolutions IR ...).